

01. DISPOSABLE CAP



Mob Cap

Material: PP
 Weight: 10- 15gsm
 Color: White, Blue, Green, Etc
 Size: 18" -24"
 Style: Single or Double Elastic
 Standard Packing: 100 pcs/ bag, 2000 pcs/ carton
 Application: Medical, Household duties, Cleaning, Beauty and Food Industry, Labs.




Bouffant Cap

Material: PP
 Weight: 10-35gsm
 Color: White, Blue, Green, Etc
 Size: 18" -24"
 Standard Packing: 100 pcs/ bag, 2000 pcs/ carton
 Application: Medical, Household duties, Cleaning, Beauty and Food Industry, Labs.



Doctor Cap

Material: PP
 Weight: 25-40gsm
 Color: Blue, Green
 Size: 13.5*64cm(length*width)
 Style: with tie on back
 Standard Packing: 100 pcs/ bag, 10 bags/ ctn, 1000 pcs/ ctn.
 Application: Industrial safety, food- processing industry, hotel, restaurant, beauty salon, clean and household disposable cover products, dental, hospital, clinic, etc



Hair Net

Material: Nylon mesh, spandex round elastic
 Weight: 2g
 Color: Black
 Size: Spandex round elastic: 4.5 cm
 Style: the edge of the cap is sewed by the serger sewing machine, so it is strong and durable
 Standard Packing: 100 pcs/ bag, 10 bag/ ctn, 1000 pcs/ ctn.
 Application: Medical supplies, industrial safety, food- processing industry, hotel, restaurant and household disposable cover products

02 FACE MASK

HAORUN MEDICAL PRODUCTS CO., Ltd



DUST MASK

-NO:BST-01V
-Material:
-200g/m2 needle-punched cotton
-Filter paper 50g/m2
-200g/m2 needle-punched cotton



DUST MASK

-NO:BST-02V
-Material:
-PP non woven 50g/m2
-Filter paper 40g/m2
-200g/m2 needle-punched cotton



DUST MASK

-NO:BST-03V
-Material:
-PP non woven 40g/m2
-Filter paper 40g/m2
-200g/m2 needle-punched cotton



Daily Use

-NO:BST-07V
-Material:
-PP non woven 40g/m2
-Active carbon non woven 50g/m2
-Filter paper 45g/m2
- Fourth filter paper 40g/m2
-PP non woven 25g/m2



Daily Use

-NO:BST-08V
-Material:
-PP non woven 40g/m2
-Filter paper 40g/m2
-Filter paper 45g/m2
- PP non woven 25g/m2



Daily Use

-NO:BST-09V
-Material:
-SMS 55g/m2
-PP non woven 40g/m2
- Filter paper FFP2 /N95 40g/m2
- PP non woven 25g/m2



DUST MASK

-NO:BST-04V
-Material:
-PP non woven 40g/m2
-Filter paper 40g/m2
-200g/m2 needle-punched cotton



DUST MASK

-NO:BST-05V
-Material:
-PP non woven 40g/m2
-Active carbon 50g/m2
-Filter paper 40g/m2
-200g/m2 needle-punched cotton



DUST MASK

-NO:BST-06V
-Material:
-PP non woven 40g/m2
-Active carbon 50g/m2
-Filter paper 40g/m2
-200g/m2 needle-punched cotton



Ear Loop Face Mask

Material:SBPP+MBPP+SBPP
Color:Blue, White ,Black, Green etc
Size:17.5*9.5cm
Style:Ear Loop
Standard Packing:50pcs/box,
40boxes/ctn, 2000pcs/ctn.
Application: Widely used for Medical/ Industrial/ Food/ Beauty/ Cleanroom/ Personal, etc.



Tie On Face Mask

Material:SBPP+MBPP+SBPP
Color:Blue, White ,Black, Green etc
Size:17.5*9.5cm
Style:Tie On
Standard Packing:50pcs/box,
40boxes/ctn, 2000pcs/ctn.
Application:Widely used for Medical/Industrial/Food/Beauty/Cleanroom/Personal, etc.



KN95

Material:
1.The first layer: PP non woven 55g/m2
2.The second layer: Melt blown Filter paper 25g/m2 (filtration efficiency ≥ 95%)
3.The third layer: Melt blown Filter paper 25g/m2 (filtration efficiency ≥ 95%)
4.The fourth layer: PP non woven 50g/m2

03 OTHER DISPOSABLE PRODUCTS

HAORUN MEDICAL PRODUCTS CO., Ltd



Machine Made Shoe Cover

Material:SBPP
Color:Blue
Size:15*40cm,17*40cm,19*40cm
Style:Anti-slip, Normal
Weight:20-40gsm
Standard Packing:10pcs/roll, 10rolls/bag, 2000pcs/ctn .
Application:Hospital, Clean rooms, labs,household, industry



Hand Made Shoe Cover

Material:SBPP
Color:Blue
Weight:20-40gsm
Size:15*40cm,17*40cm,19*40cm
Style:Anti-slip, Normal
Standard Packing:10pcs/roll, 10rolls/bag, 2000pcs/ctn .
Application:Hospital, Clean rooms, labs, household, industry



CPE Shoe Cover

Material:CPE
Color:Blue
Weight:2-5g/pcs
Size:16*38cm,17*42cm
Type:Machine made
Standard Packing:10pcs/roll, 10rolls/bag, 2000pcs/ctn .
Application:Hospital, Clean rooms, labs, household, industry



Beard Cover

Material:SBPP
Color:White
Weight:10-20gsm
Size:10-24"
Style:Machine made
Standard Packing:100pcs/bag, 2000pcs/ctn .
Application:Food Industry, Clean rooms, labs, household, industry



Over Sleeves

Material:PE
Color:White,Blue
Weight:2-5g/pcs
Size:22*46cm
Style:Machine made
Standard Packing:100pcs/bag, 2000pcs/ctn .
Application:Food Industry, Clean rooms, labs, household, industry



CPE Gown

Material:CPE
Color:White,Blue
Weight:30-50g/pcs
Size:S-XXL
Style:Long sleeves,opening cuff with thumb loop,tie on waist .
Standard Packing:1pcs/bag, 10pcs/big bag, 100pcs/ctn .
Application:Food Industry, Clean rooms, labs, household, industry

BED SHEET/ BED COVER



Material:PP/SMS
Color:White,Blue
Weight:20-45g/pcs
Size:80*180cm, 140*220cm,120*220cm
Standard Packing:10 pcs per bag, 10 bags per carton .
Application:home textile,medical use,industry

NON- WOVEN FABRIC

Eco-friendly,water repellent
Anti- UV(1%-5%), anti-bacteria,anti-static,
flame retardant function as request
Tear resistant, shrink-resistant
Strong strength and elongation,soft, non-toxic
Excellent property of air permeability



04. DISPOSABLE CLOTHES

Material: PP/ SMS/ Microporous
Color: White/ Blue/ Red/ Orange/ Grey, etc
Weight: 30-65 gsm
Garment features: comfortable, dust proof, light weight
Size: S-3XL
Sewing: Serged Seam
Package: 1 pcs/ bag, 50 pcs/ Carton
Application: General purpose, Clean room, Construction, Industry, Mining, Agriculture, Asbestos Stripping, Non-hazardous Operations, Pharmaceutical

Style:
Conjoined, or seam collection
With or without hood
With or without boots
With or without elastic on the waist
Elastic cuff or knitted cuff
With or without color tape on the flap
With or without flap at front
With or without ring
One way or two way zipper



With Seal tape
microporous coverall



PP Coverall



SMS Coverall



Microporous Coverall



SURGICAL GOWN

- Material: 30~65gsm SMS/SMMS
- Size: S /M / L /XL /2XL
- Style: with knitted cuff, tie on neck or velcro on neck, 2 waist belts or 4 waist belts.
With/ without Tissue paper/ wrap cloth
- process: A. Sew B. Ultrasonic
- Package: 1 pc/sterile pouch, 50 pcs/carton



ISOLATION GOWN

- Material: 20~50gsm PP/ SMS/PP+PE
- Size: S /M / L /XL /2XL
- Color: blue, white, yellow, green
- Style: elastic cuff/knitted cuff, tie on neck and waist. Non sterile.
- Package: 10 pcs/bag, 10 bags/carton



LAB COAT/WORKING SUIT

- Material: 20~50gsm PP / SMS/PE
- Size: S /M / L /XL /2XL /3XL
- Color: white, blue, red.
- Style: with elastic/knitted cuff, velcro/snaps/ zipper in front , T-shirt/Knitted/Korea collar, with or without pockets.
- Package: 1 pc/ bag, 100 pcs/ carton



Scrub Suits

- Material: 30~65gsm PP/ SMS/PP+PE
- Size: S /M / L /XL /2XL /3XL
- Color: blue, green
- Style: V/round color,short sleeves or non sleeves, with or without pockets.
- Package: 1 pcs/bag, 50 bags/carton