



普通型气管插管

Endotracheal Tube



说明: FEATURES:

- 无毒医用级PVC制造, 透明、柔软、光滑。
- 大容量低压力气囊。
- X光造影线贯穿管身。
- 带墨菲孔。
- Made of non-toxic medical-grade PVC, transparent, soft and smooth.
- With high volume low press cuff.
- Radio opaque line through the length for x-ray visualization.
- Available with Murphy eye only.

产品编号 Product No.	规格 (#) Size (#)	产品编号 Product No.	规格 (#) Size (#)	备注 Remark
ET20PC	2.0	ET65PC	6.5	分带囊 和无囊 cuffed and uncuffed
ET25PC	2.5	ET70PC	7.0	
ET30PC	3.0	ET75PC	7.5	
ET35PC	3.5	ET80PC	8.0	
ET40PC	4.0	ET85PC	8.5	
ET45PC	4.5	ET90PC	9.0	
ET50PC	5.0	ET95PC	9.5	
ET55PC	5.5	ET100PC	10.0	
ET60PC	6.0			

加强型气管插管

Reinforced Endotracheal Tube



说明: FEATURES:

- 无毒医用级PVC制造, 透明、柔软、光滑。
- 大容量低压力气囊。
- 通气管经过弹簧加强, 可防止压扁、扭曲。
- 适合病人各种体位, 特别卧位手术。
- Made of non-toxic medical-grade PVC, transparent, soft and smooth.
- With high volume low press cuff.
- Spiral reinforcement minimizes crushing or kinking.
- Conform to any patient posture, especially to operation of decubitus.

产品编号 Product No.	规格 (#) Size (#)	产品编号 Product No.	规格 (#) Size (#)	备注 Remark
ET30PCR	3.0	ET65PCR	6.5	分带囊 和无囊 cuffed and uncuffed
ET35PCR	3.5	ET70PCR	7.0	
ET40PCR	4.0	ET75PCR	7.5	
ET45PCR	4.5	ET80PCR	8.0	
ET50PCR	5.0	ET85PCR	8.5	
ET55PCR	5.5	ET90PCR	9.0	
ET60PCR	6.0	ET95PCR	9.5	

双腔硅胶导尿管

2-Way Silicone Foley Catheter



说明: FEATURES:

- 采用进口医用硅胶制作。
- 球囊具有良好的均衡性和拉伸性, 使用过程中更加安全、可靠。
- 长度: 330mm或400mm。
- Made of medical-grade silicone.
- Symmetrical balloon expands equally in all directions so as to perform its function safely and efficiently.
- Length: 330mm or 400mm.

产品编号 Product No.	规格 (Fr) Size (Fr)	色标号 Color Code	球囊容量(mL) Balloon Volume(mL)	备注 Remark
FC0615S2	6	粉红色/Pink	3	有导丝 With stylet
FC083/5S2	8	浅蓝色/Light Blue	3-5	
FC103/5S2	10	黑色/Black	3-5	
FC1215/30S2	12	白色/White	5-15	无导丝 Without stylet
FC1415/30S2	14	绿色/Green	5-15/30	
FC1615/30S2	16	橙色/Orange	5-15/30	
FC1815/30S2	18	红色/Red	5-15/30	
FC2015/30S2	20	黄色/Yellow	5-15/30	
FC2215/30S2	22	紫色/Purple	5-15/30	
FC2415/30S2	24	蓝色/Blue	5-15/30	
FC2615/30S2	26	粉红色/Pink	5-15/30	

三腔硅胶导尿管

3-Way Standard Silicone Foley Catheter



说明: FEATURES:

- 采用进口医用硅胶制作。
- 球囊具有良好的均衡性和拉伸性, 使用过程中更加安全、可靠。
- 长度: 400mm。
- Made of medical-grade silicone.
- Symmetrical balloon expands equally in all directions so as to perform its function safely and efficiently.
- Length: 400mm.

产品编号 Product No.	规格 (Fr) Size (Fr)	色标号 Color Code	球囊容量(mL) Balloon Volume(mL)	备注 Remark
FC1630S3	16	橙色/Orange	30	3腔 3-way
FC1830S3	18	红色/Red	30	
FC2030S3	20	黄色/Yellow	30	
FC2230S3	22	紫色/Purple	30	
FC2430S3	24	蓝色/Blue	30	





硅胶喉罩

Silicone Laryngeal Mask

说明: FEATURES:

- 采用进口医用级硅胶原料、无毒无刺激。
- 罩囊部选用柔软的医用级硅胶材料, 外形设计与咽喉部轮廓一致, 最大限度减少对患者的刺激, 提高了与喉部的密封性。
- 可适用于成人、儿童、新生儿。
- Made of medical-grade silicone.
- The cuff is made of medical-grade silicone, and its specially designed shape coincides with the laryngopharynx well, reducing stimulation to patient body and improving the cuff seal.
- Suitable for adult, children and infant use.

产品编号 Product No.	规格 (#) Size (#)	色标号 Color Code	病人体重(KG) Patient Weight(KG)	备注 Remark
LM10S1	1.0	黑色/Black	0-5	一次性使用 Disposable
LM15S1	1.5	白色/White	5-10	
LM20S1	2.0	绿色/Green	10-20	
LM25S1	2.5	橙色/Orange	20-30	
LM30S1	3.0	红色/Red	30-50	
LM40S1	4.0	黄色/Yellow	50-70	
LM50S1	5.0	紫色/Purple	70-100	



PVC喉罩

PVC Laryngeal Mask

说明: FEATURES:

- 通气管采用医用级PVC原料。
- 罩囊部选用柔软的医用级硅胶材料, 外形设计与咽喉部轮廓一致, 最大限度减少对患者的刺激, 提高了与喉部的密封性。
- 可适用于成人、儿童、新生儿。
- The tube is made of medical-grade PVC.
- The cuff is made of medical-grade silicone, and its specially designed shape coincides with the laryngopharynx well, reducing stimulation to patient body and improving the cuff seal.
- Suitable for adult, children and infant use.

产品编号 Product No.	规格 (#) Size (#)	色标号 Color Code	病人体重(KG) Patient Weight(KG)	备注 Remark
LM10P1	1.0	黑色/Black	0-5	一次性使用 Disposable
LM15P1	1.5	白色/White	5-10	
LM20P1	2.0	绿色/Green	10-20	
LM25P1	2.5	橙色/Orange	20-30	
LM30P1	3.0	红色/Red	30-50	
LM40P1	4.0	黄色/Yellow	50-70	
LM50P1	5.0	紫色/Purple	70-100	



加强型喉罩

Reinforced Laryngeal Mask

说明: FEATURES:

- 通气管经钢丝加强, 可防止扭结、或压扁。
- 管体光滑柔软、透明, 抗折曲性极佳。
- 可适用于成人、儿童、新生儿。
- Spiral reinforcement minimizes crushing or kinking.
- Smooth, transparent and kink-resistant tube.
- Suitable for adult, children and infant use.

产品编号 Product No.	规格 (#) Size (#)	色标号 Color Code	病人体重(KG) Patient Weight(KG)	备注 Remark
LM10S1R	1.0	黑色/Black	0-5	一次性使用 Disposable
LM15S1R	1.5	白色/White	5-10	
LM20S1R	2.0	绿色/Green	10-20	
LM25S1R	2.5	橙色/Orange	20-30	
LM30S1R	3.0	红色/Red	30-50	
LM40S1R	4.0	黄色/Yellow	50-70	
LM50S1R	5.0	紫色/Purple	70-100	



可重复性使用硅胶喉罩

Reusable Silicone Laryngeal Mask

说明: FEATURES:

- 采用进口医用级硅胶原料, 无毒无刺激。
- 可适用于成人、儿童、新生儿。
- 高压蒸汽灭菌, 可重复使用40次。
- Made of medical-grade silicone.
- Suitable for adult, children and infant use.
- Autoclave sterilization only; Can be reused up to 40 times.

产品编号 Product No.	规格 (#) Size (#)	色标号 Color Code	病人体重(KG) Patient Weight(KG)	备注 Remark
LM10S2	1.0	黑色/Black	0-5	可重复性使用 Reusable
LM15S2	1.5	白色/White	5-10	
LM20S2	2.0	绿色/Green	10-20	
LM25S2	2.5	橙色/Orange	20-30	
LM30S2	3.0	红色/Red	30-50	
LM40S2	4.0	黄色/Yellow	50-70	
LM50S2	5.0	紫色/Purple	70-100	





吸痰管

Suction Tube



说明: FEATURES:

- 采用医用级PVC原料, 柔软透明。
- 不同型号有颜色区分。
- Made of non-toxic medical-grade PVC, clear and soft.
- Color-coded connector for identification of different sizes.

产品编号 Product No.	规格 (Fr) Size (Fr)	色标号 Color Code	长度(mm) Length(mm)
ST5038S	5	灰色/Gray	380
ST6038S	6	浅绿色/Light Green	380
ST8050S	8	浅蓝色/light blue	500
ST1050S	10	黑色/Black	500
ST1250S	12	白色/White	500
ST1450S	14	绿色/Green	500
ST1650S	16	橙色/Orange	500
ST1850S	18	红色/Red	500

口咽通气道

Guedel Airway



说明: FEATURES:

- 半硬无毒聚乙烯材料制成, 带色标。
- 中央单通道。
- Made from semirigid non-toxic PE color bide back.
- Guedel style enclosed center channel.

产品编号 Product No.	规格 (#) Size (#)	色标号 Color Code	长度(mm) Length(mm)
OA40S1	4.0	粉红色/Pink	40
OA50S1	5.0	蓝色/Blue	50
OA60S1	6.0	黑色/Black	60
OA70S1	7.0	白色/White	70
OA80S1	8.0	绿色/Green	80
OA90S1	9.0	黄色/Yellow	90
OA100S1	10.0	红色/Red	100
OA110S1	11.0	浅蓝色/light blue	110
OA120S1	12.0	橙色/Orange	120

胃管

Stomach Tube



说明: FEATURES:

- 医用PVC或硅胶制造。
- X光造影线贯穿管身。
- Made from medical grade PVC or Silicone.
- Radio-opaque line through the tube.

产品编号 Product No.	规格 (Fr) Size (Fr)	色标号 Color Code	长度(mm) Length(mm)
ST05P	5	灰色/Gray	50
ST06P	6	浅绿色/Light Green	110
ST08P	8	浅蓝色/light blue	110
ST10P	10	黑色/Black	130
ST12P	12	白色/White	130
ST14P	14	绿色/Green	130
ST16P	16	橙色/Orange	130
ST18P	18	红色/Red	130

氧气面罩

Oxygen Mask



说明: FEATURES:

- 透明、柔软PVC制成。
- 可调鼻夹。
- 可配置210cm的供氧管。
- Made of clear soft PVC.
- Adjustable nose clip assures comfortable fit.
- 210CM Oxygen supply tube can be configured.

产品编号 Product No.	型号 Type
102MM	成人带管/Adult, with tubing
103MM	成人不带管/Adult, without tubing
102MM-E	成人加长带管/Adult, elongated with tubing
103MM-E	成人加长不带管/Adult, elongated without tubing
222MM	小儿, 带管/Pediatric, with tubing
223MM	小儿, 不带管/Pediatric, without tubing
222MM-E	小儿加长带管/Pediatric elongated, with tubing
223MM-E	小儿加长不带管/Pediatric elongated, without tubing





带抽吸气管插管(普通型) Endotracheal Tube with Aspiration Lumen

说明: FEATURES:

- 可有效抽吸聚集在声门下分泌物;
- 抽吸分泌物时不会对声道粘膜造成损伤;
- 抽吸方式便捷安全;
- 有助于降低医院内呼吸机相关性肺炎的发生率;
- 延缓呼吸机相关性肺炎的早期发生。
- Effective suctioning of accumulated secretions in the subglottic space;
- No trauma to vocal cords when suctioning;
- Convenient and safe;
- Helping to reduce the rate of VAP(Ventilator-Associated Pneumonia)in the hospital;
- Helping delay the early onset of VAP.

规格/Size (I.D.mm)	圆形头部,有囊 Rounded beveled tip,with Cuff	斜面头部,有囊 Beveled tip,with Cuff	包装/Package
5.0	202050	202150	10 pcs/box,10 boxes/ctn
5.5	202055	202155	10 pcs/box,10 boxes/ctn
6.0	202060	202160	10 pcs/box,10 boxes/ctn
6.5	202065	202165	10 pcs/box,10 boxes/ctn
7.0	202070	202170	10 pcs/box,10 boxes/ctn
7.5	202075	202175	10 pcs/box,10 boxes/ctn
8.0	202080	202180	10 pcs/box,10 boxes/ctn
8.5	202085	202185	10 pcs/box,10 boxes/ctn
9.0	202090	202190	10 pcs/box,10 boxes/ctn



带抽吸气管插管(加强型) Endotracheal Tube with Aspiration Lumen

说明: FEATURES:

- 可有效抽吸聚集在声门下分泌物;
- 抽吸分泌物时不会对声道粘膜造成损伤;
- 抽吸方式便捷安全;
- 有助于降低医院内呼吸机相关性肺炎的发生率;
- 延缓呼吸机相关性肺炎的早期发生。
- Effective suctioning of accumulated secretions in the subglottic space;
- No trauma to vocal cords when suctioning;
- Convenient and safe;
- Helping to reduce the rate of VAP(Ventilator-Associated Pneumonia)in the hospital;
- Helping delay the early onset of VAP.

规格/Size (I.D.mm)	圆形头部,有囊 Rounded beveled tip,with Cuff	斜面头部,有囊 Beveled tip,with Cuff	包装/Package
5.0	202050	202150	10 pcs/box,10 boxes/ctn
5.5	202055	202155	10 pcs/box,10 boxes/ctn
6.0	202060	202160	10 pcs/box,10 boxes/ctn
6.5	202065	202165	10 pcs/box,10 boxes/ctn
7.0	202070	202170	10 pcs/box,10 boxes/ctn
7.5	202075	202175	10 pcs/box,10 boxes/ctn
8.0	202080	202180	10 pcs/box,10 boxes/ctn
8.5	202085	202185	10 pcs/box,10 boxes/ctn
9.0	202090	202190	10 pcs/box,10 boxes/ctn



PVC标准型导尿管 PVC Nelaton Catheter

说明: FEATURES:

- 医用PVC制造
- 长度:400mm
- Made from medical grade
- Length:400mm

产品编号 Product No.	规格 (#) Size (#)	色标 Color Code
121106	6	淡绿色/Light Green
121108	8	淡蓝色/Light Blue
121110	10	黑色/Black
121112	12	白色/White
121114	14	绿色/Green
121116	16	橙色/Orange
121118	18	红色/Red
121120	20	黄色/Yellow



双腔支气管插管 Double lumen endobronchial intubation

说明: FEATURES:

- 示压气阀系列可方便指示囊内压力,对病人气道损伤小
- 低压高容气囊
- 进口PVC料
- 纯硅胶内芯单向阀可用于核磁共振检查
- Pressure gauge valve series can be convenient for indicating the intracapsular pressure, airway injury to small
- Low pressure high volume set of sac
- Imports of PVC material
- Pure silica core one-way valve can be used in magnetic resonance imaging

型号	规格	左侧	右侧	
			S型囊	双囊
普通型	柱状囊	26#、28#、31#、33#、35#、37#、39#、41#	33#、35#、37#、39#、41#	26#、28#、31#、33#、35#、37#、39#、41#
	橄榄状球囊	26#、28#、31#、33#、35#、37#、39#、41#	33#、35#、37#、39#、41#	26#、28#、31#、33#、35#、37#、39#、41#
接头固定	柱状囊	35#、37#、39#、41#	35#、37#、39#、41#	35#、37#、39#、41#
	橄榄状球囊	35#、37#、39#、41#	35#、37#、39#、41#	35#、37#、39#、41#
带隆突钩型	柱状囊	35#、37#、39#、41#	35#、37#、39#、41#	35#、37#、39#、41#
	橄榄状球囊	35#、37#、39#、41#	35#、37#、39#、41#	35#、37#、39#、41#






双腔乳胶导尿管 2-Way Silicone Foley Catheter

说明: FEATURES:

- 采用天然乳胶制作, 100%硅胶涂层。
- 球囊具有良好的均衡性和拉伸性, 使用过程中更加安全、可靠。
- 长度: 330mm或400mm。
- Made of latex, 100% silicone coated.
- Symmetrical balloon expands equally in all directions so as to perform its function safely and efficiently.
- Length: 330mm or 400mm.

产品编号 Product No.	规格 (Fr) Size (Fr)	色标号 Color Code	球囊容量(mL) Balloon Volume(mL)	备注 Remark
LC0615S2	6	粉红色/Pink	3	有导丝 With stylet
LC083/5S2	8	浅蓝色/Light Blue	3-5	
LC103/5S2	10	黑色/Black	3-5	
LC1215/30S2	12	白色/White	5-15	
LC1415/30S2	14	绿色/Green	5-15/30	无导丝 Without stylet
LC1615/30S2	16	橙色/Orange	5-15/30	
LC1815/30S2	18	红色/Red	5-15/30	
LC2015/30S2	20	黄色/Yellow	5-15/30	
LC2215/30S2	22	紫色/Purple	5-15/30	
LC2415/30S2	24	蓝色/Blue	5-15/30	
LC2615/30S2	26	粉红色/Pink	5-15/30	



三腔硅胶导尿管 3-Way Latex Silicone Foley Catheter

说明: FEATURES:

- 采用天然乳胶制作, 100%硅胶涂层。
- 球囊具有良好的均衡性和拉伸性, 使用过程中更加安全、可靠。
- 长度: 400mm。
- Made of latex, 100% silicone coated.
- Symmetrical balloon expands equally in all directions so as to perform its function safely and efficiently.
- Length: 400mm.

产品编号 Product No.	规格 (Fr) Size (Fr)	色标号 Color Code	球囊容量(mL) Balloon Volume(mL)	备注 Remark
LC1630S3	16	橙色/Orange	30	3腔 3-way
LC1830S3	18	红色/Red	30	
LC2030S3	20	黄色/Yellow	30	
LC2230S3	22	紫色/Purple	30	
LC2430S3	24	蓝色/Blue	30	



扬克头吸引管 Yankauer Handle And Suction Connecting Tube

说明: FEATURES:

- K料或亚克力制作、PVC连接管。
- 规格齐全, 大小头, 有孔无孔均有。
- Made of K-resin and medical PVC.
- Supporting the use.

产品编号 Product No.	规格 (Fr) Size (Fr)	颜色 Color Code	备注 Remark
SC001A	1/4"-1.8m	Light Blue	
SC001B	1/4"-1.8m	Light Blue	
SC001A	1/16"-1.8m	Light Blue	
SC001B	3/16"-3.8m	Light Blue	
YH001A	Crown to without vent	Light Blue	
YH001B	Crown to with vent	Light Blue	
YH002A	Crown to without vent	Light Blue	
YH002B	Standnrd with vent	Light Blue	



鼻咽通气道 Nasopharyngeal airway

说明: FEATURES:

- PVC制作, 柔软舒适。
- 仅供鼻腔通气用。
- 尾端有喇叭口型、平口型、平口加厚型和滑动环型供选择。

产品编号 Product No.	规格 (#) Size (#)	产品编号 Product No.	规格 (#) Size (#)
NA20PC	2.0	NA25PC	2.5
NA30PC	3.0	NA35PC	3.5
NA40PC	4.0	NA45PC	4.5
NA50PC	5.0	NA55PC	5.5
NA60PC	6.0	NA65PC	6.5
NA70PC	7.0	NA75PC	7.5
NA80PC	8.0	NA85PC	8.5
NA90PC	9.0		





气管插管固定器 Endotracheal Tube Holder

说明: FEATURES:

- 适用于各种气管插管。
- 产品仅供一次性使用。
- Fits various sizes of endotracheal tubes, from Id5.5 to Id10.
- For single use only.

产品编号 Product No.	规格 (Fr) Size (Fr)	色标号 Color Code	长度(mm) Length(mm)
EF1E01	A		380
EF1E02	B		380



鼻氧管 Nasal Oxygen Cannula

说明: FEATURES:

- 由无毒PVC制成。
- 包装规格: 大、中、小。
- 包装材料: 塑料袋或纸塑料袋包装。
- 产品仅供一次性使用。
- Made from non-toxic PVC.
- Packaging
- Available with single-hole nare and twin-hole nare.
- Anti-crush oxygen tubing.

产品编号 Product No.	规格 Size	描述说明 DESCRIPTION	备注 Remark
C0201	Adult	standard tip	一次性使用 Disposable
C0202	Pediatric	standard tip	
C0203	Infant	standard tip	
C0204	Adult	flared tip	
C0205	Adult	soft tip	
C0206	Pediatric	soft tip	
C0207	Infant	soft tip	
C0208	Neonate	soft tip	
C0209	Adult	curved tip	



导丝 Intubating Stylet

说明: FEATURES:

- 医用PVC或硅胶制造。
- Made from medical grade PVC or Silicone.

产品编号 Product No.	规格 (Fr) Size (Fr)	长度(mm) Length(mm)	备注 Remark
Is061	6	305	
Is101	10	380	
Is141	14	380	



气管切开插管 Tracheostomy Tube

说明: FEATURES:

- 由医用级PVC制成。
- 包装材料: 泡壳包装。
- 产品仅供一次性使用。
- Made from medical grade PVC.
- Supplied sterile in individual peelable blister pack.
- For single use only.

产品编号 Product No.	规格 (#) Size (#)	产品编号 Product No.	规格 (#) Size (#)	备注 Remark
ETT30PC	3.0	ETT65PC	6.5	分带囊 和无囊 cuffed and uncuffed
ETT35PC	3.5	ETT70PC	7.0	
ETT40PC	4.0	ETT75PC	7.5	
ETT45PC	4.5	ETT80PC	8.0	
ETT50PC	5.0	ETT85PC	8.5	
ETT55PC	5.5	ETT90PC	9.0	
ETT60PC	6.0	ETT95PC	9.5	





全麻气管插管包 Anaesthesia Windpipe Intubation Package

说明: FEATURES:

- 适用于各种气管插管。
- 产品仅供一次性使用。
- Fits various sizes of endotracheal tubes, from Id5.5 to Id10.
- For single use only."

产品编号 Product No.	规格 Size
AWIP01	配加强带囊插管ID2.0-10.0 Equipped with reinforced endotracheal tube with cuffed ID2.0-10.06
AWIP02	配加强无囊插管ID2.0-10.0 Equipped with reinforced endotracheal tube without cuffed ID2.0-10.0
AWIP03	配普通带囊插管ID2.0-10.0 Equipped with endotracheal tube with cuffed ID2.0-10.0
AWIP04	配普通无囊插管ID2.0-10.0 Equipped with endotracheal tube without cuffed ID2.0-10.0



麻醉面罩 Anesflex Mask

说明: FEATURES:

- 麻醉面罩根据面部结构设计完美。
- 密封垫块柔软提供最大的舒适度。
- 无毒,无害,强度高。
- 波浪外形的缓冲。
- 7种可用大小,以适应各种病人。所有尺寸均提供一个定制适合性和舒适性。
- 22毫米彩色指环利于手持移动。
- 麻醉面罩可用于麻醉、呼吸机、氧气呼吸机、高压室、无痛分娩机以及救援呼吸袋。
- 适合使用于麻醉、呼吸系统和人工呼吸器。

- Anesflex mask is designed to fit perfectly for the face structure.
- Airtight cushion pat with superior softness provides maximum comfort.
- Non-toxic, non-odor, high strength.
- Contoured cushion creates a gentle with the patient's face.
- Available in 7 sizes to accommodate a full range of patients. All sizes include an inflation port providing a customized fit and comfort.
- 22mm color-coded retaining rings are removable for hand-held procedures.
- Anesflex mask can be connected with the anesthesia, breathing machine, oxygen breathing machine, high pressure chamber, imbibe-way painless childbearing machine as well as the rescue respiration bag.
- Suitable for use of anesthesia, respiratory and resuscitator.

产品编号 Product No.	规格 (#) Size (#)	产品编号 Product No.	规格 (#) Size (#)
Am00	0-Pre neonate	Am04	4-Small Adult
AM01	1-Neonate	AM05	5-Medium Adult
AM02	2-Infant	AM06	6-Large Adult
AM03	3-Child		



HAORUN MEDICAL PRODUCTS CO.,LTD

TEL 0086-574-87706193 WHATSAPP/WECHAT/TEL:008613566542343
EMAIL:SALES@HAORUNMED.COM WWW.HAORUNMED.COM